C TECNIC



ePilot

bioreactor single use

ePilot tangential flow filtration

Multi use

Vessels

Single use Vessels

Cassette: 0.5-5 Hollow Fiber: 0.1-6.5

Filter area (m²)

1200x2020x910

Size in mm

600x1687x595 Size in mm

1546x1412x524

Size in mm

 $10 \cdot 20 \cdot 30 \cdot 50$ Volumes (L)

30 . 50 Volumes (L) 50 - 200

Volumes (L)

Stainless steel

Housing Material

Stainless steel

Housing Material

Versatile Holder Accommodating

Range of filter options

1x Variable Speed

eSCADA Advanced

Software

eSCADA Advanced

1x Fixed Speed Peristaltic Pump

4x Variable Speed

Servomotor 400-750

Motor (W)

Peristaltic Pumps

4x Integrated Variable Pumps 3x External Variable Pumps

Software

Peristaltic Pumps

Compatibility

TECNIC eBAG 3D STR

20 LPM Recirculation Pump

230VAC | 50Hz | 9A

Power Supply

ePilot Pilot Scale Solutions

Optimizing bioprocessing outcomes with innovative engineering and precise regulation for pilot plants.



C TECNIC

STR ePilot BR

Microbial / Cellular · 10 L up to 50 L

High-performance automatic bioreactor with cellular and microbial configurations, an excellent choice for pilot scale development and production.



STR Single Use ePilot STR SU

Dual Purpose · 30 L and 50 L

A scalable, dual-purpose bioreactor featuring precise temperature control, efficient agitation and wide instrumentation, designed for TECNIC's single use eBAG.



TEE ePilot TFF

Hollow Fiber / Cassette · Up to 6.5 m² of surface area

Automatic pilot Tangential Flow Filtration system with high efficiency and versatility for bioprocessing, featuring options for concentration, diafiltration, and membrane surface adjustments.

Global reach, local support

Our global distributor network is ready to provide specialized products and support in your region. Please contact your local distributor to work with TECNIC

sales@tecnic.eu

















Adaptable and interconnected

Build a pilot plant that grows with your needs. ePilot and ePlus equipment work hand-in-hand to offer scalability, efficiency, and flexibility. Transition from simple processes to advanced setups without compromising process control or quality.

ePilot bioreactor

Scalable pilot bioreactor for process development with precise control over key parameters.

ePilot bioreactor single use

Bioreactor with single-use bags for flexible, contamination-free workflows.

ePilot tangential flow filtration

A tangential flow filtration system for pilot-scale operations, perfect for concentration and diafiltration processes. Designed for ease of use and high-performance membrane filtration.





